






Olympus® SZX7 System Diagram



*1 Focusing unit (SZX2-FOF, SZX-FOFH or SZX-FO) and SZX-STF are required when mounting fluorescent unit.
 *2 SZX2-FOFH and SZX2-STL2 are required when using SZX-SDO2.
 *3 Please contact your nearest Olympus dealer for applicable cameras.
 *4 SZH-P400 and SZH-P600 can be attached to the transmitted light illuminators.


*5 Equipped to SZX2-ILLC10.
 *6 Please contact your nearest Olympus dealer for applicable observation tube combination.
 *7 Not available in some areas.
 *8 SZ-LW61 cannot be combined with SZX2-ILST.

Olympus® SZX61/51 System Diagram




*1 SZH-P400 and SZH-P600 can be attached to the transmitted light illuminators. *2 Not available in some areas.


Compact fiber optics system/SZ2-LGB



High power fiber optics system/LG-PS2




Ring light illumination



*3 SZ-LW61 cannot be combined with SZ2-ILST.

Reflected illumination/LSGA



*4 IHA is mountable to SZ2-ST, SZ2-STB1, SZ2-STB3, SZ2-ST5 and TL-3.

*5 IHD is mountable to SZ2-STB1 and SZ2-ST5.

*6 IHE is mountable to SZ2-ST.